



## LIGHTFORCE 2.0 FEATURES

# Custom Tubes

### **LightForce 2.0 Feature: Custom Tubes**

#### **What is it?**

LightForce 2.0 introduces Custom Tubes for 6s and 7s.

#### **What problem does it solve?**

These new custom tubes function like all other LightForce brackets and truly makes the LightForce system fully customizable.

#### **How do I get it?**

Custom tubes on 6s will be available for new cases submitted after April 7. Custom tubes on 7s will be available over the next few weeks after the launch of LF 2.0. You will be notified once you have access to custom tubes on 7s.

#### **How does it work?**

Custom tubes will be planned for all cases where there is adequate anatomy and scan data to support custom tubes. The scan should include all tooth anatomy and a minimum of 0.5 mm of gingiva surrounding the terminal molar. The first molar bracket will be included in the posterior IDB tray segment, same as before, and the second molar brackets will be provided in a single tooth tray.

#### **How does it impact my practice?**

Inclusion of 6s and 7s allow you to create a more comprehensive treatment plan with increased vertical and angular control of molars. Custom tubes have better adaptation to molar anatomy for increased bond strength. Custom tubes also provide the added benefit of a fully delegated bonding and rebonding process.

#### **Can I use the “second scan” option to obtain custom tubes for 6s and 7s?**

For cases created after access to 7s in early May, we will allow you to request custom tubes for 6s and 7s that were not fully erupted in the initial scan. We cannot guarantee that we will be able to fulfill this request, as we are constrained by the bracket placement of the original LightPlan, but we will make a best effort to support this. **Custom tubes for 6s and 7s will not be available for cases created prior to full release in early May.**



## LIGHTFORCE 2.0 FEATURES

# Custom Tubes

### Additional FAQs

**Q: For cases that are currently in Planning, can they add custom 6s?**

A: If the doctor is not part of the beta, custom 6s are not available for any cases created before April 7. Cases submitted from that date on will have custom 6s.

**Q: If Dr does not take a scan of 7's and does not want 7's included will the case still be put on hold?**

A: No, the case should not be put on hold if the doctor indicates that they don't want 7s in the Rx form.

**Q: When are custom tubes available?**

A: After April 7, Custom 6s will be available for all practices for cases created after that date. Practices will be notified on a rolling basis through early May when 7s are available to them, when they are expected to be GA.

**Q: For cases submitted prior to when I have access to 7s, can I use mixed dentition to get 7s later on?**

A: No, mixed dentition will only be available for 7s for cases created after 7s are released in early May. They will not be available for historical cases via Mixed Dentition/Second Scans.